











<b>Type A</b> <u>US Standard</u> Beckman, Corning, Orion, Fisher etc.	<b>Type B</b> <u>BNC</u> IQ, Chemtrix, Cole-Parmer, Consort, Hach, Beckman, Jenco, Hanna, Corning, recent Orion	<b>Type C</b> <u>Radiometer</u> (after '74/'75) Crison	<b>Type D</b> <u>Pin</u> Sargent-Welch models EP & RB, Fisher 150, most US reference electrodes	<b>Type E</b> <u>Photovolt</u>	<b>Type F</b> <u>Miniature Phone Plug A</u> Corning Models 5&6
					
<b>Type G</b> <u>Miniature Phone Plug B</u> Corning Model 3	<b>Type H</b> <u>Spade Lug</u> Industrial Electrodes	<b>Type I</b> <u>DIN 19262</u> German meters Metrohm type E, Mettler, Consort & recent Polymetron	<b>Type K</b> <u>Belling Lee</u> British meters EIL, PYE-Unicam	<b>Type L</b> <u>Beckman</u> Beckman Chem-Mate, Pistol, Expand-mate & Select-Mate	<b>Type M</b> <u>Electrofact</u> Electrofact & some Philips
					
<b>Type N</b> <u>Lemo</u> Seybold	<b>Type R</b> <u>Metrohm</u> recent Metrohm	<b>Type T</b> <u>No Connector</u> Stripped & tinned	<b>Type U</b> <u>Banana</u> Reference & Redox electrodes	<b>Type V</b> <u>Lemo</u> recent Mettler	





pHoenix Electrode Sales CC ! P.O. Box 5487 ! Cresta ! 2118  
Tel. (011) 792 1210 ! Fax. (011) 793 6863 ! e-mail. sales@ph.co.za