


Conductivity Electrodes - Laboratory Immersion Type



These conductivity electrodes are designed to be used with portable or laboratory instruments. Platinum cells are platinized platinum.

When ordering specify temperature compensator, if any, cable length, and connector type compatible with your conductivity meter.

Part No.	Body	Electrode	Temperature	Cell Constants (cm-1)
5731301	Epoxy	Graphite	-5 to +80°C	0.1, 0.5, 1.0, 2.0, 5.0, 10.0
5771302	Glass	Platinum 2 bands	-5 to +110°C	0.1, 1.0, 2.0
5771303	Glass	Platinum 3 bands	-5 to +110°C	0.1, 0.5, 1.0
5771304	Glass	Platinum 2 bands	-5 to +110°C	0.1, 0.5, 1.0, 2.0
5771305	Glass	Platinum 3 bands	-5 to +110°C	0.1, 0.5, 1.0, 2.0



pHoenix Electrode Sales CC ! P.O. Box 5487 ! Cresta ! 2118
 Tel. (011) 792 1210 ! Fax. (011) 793 6863 ! e-mail. sales@ph.co.za