

Glass Body Laboratory Half Cells

Although the modern trend is towards the use of combination electrodes for pH, ORP and ion selective measurements there is still a requirement for half cells.


pHoenix Electrode Company reference and measuring electrodes give excellent performance and are economically priced.



5773301	Glass Mono pH	0~13pH	GX-2	0~100°C	
	Water, macro samples, TRIS buffer, sea water				12mm dia. x 110mm



5771423	Glass Ag/AgCl reference	Refillable		0~100°C	Single Junction Ag/AgCl
	Water, macro samples				12mm dia. x 110mm



5771424	Glass Calomel Reference	Refillable		-5~80°C	Single Junction Calomel
	Lab flasks, test tubes				12mm dia. x 110mm



pHoenix Electrode Sales CC ! P.O. Box 5487 ! Cresta ! 2118
 Tel. (011) 792 1210 ! Fax. (011) 793 6863 ! e-mail. sales@ph.co.za